

# ROLLATOR NEW DESIGN





## Rollator New Design - W2230 002 001

#### **Product description**

This stylish, lightweight, compact rollator has been studied to provide comfortable support during everyday outings. Easy to handle, stable and light, the rollator New Design provides a natural walking position. Let's take a safe and comfortable break thanks to the parking brakes, the cushioned seat and flexible backrest to lean on without falling over backwards. Very convenient, it is equipped with a removable bag to go shopping or to store personal belongings. Totally foldable, it is ultra-compact and fits easily into the boot of a small car.

This product is CE marked according to the European council directive 93/42/EEC and its relatives or the Regulation 2017/745 on medical devices.








Ultra-compact when fully folded.

### **General specifications**


- Foldable rollator in aluminum.
- Removable shopping bag and left/right reversible cane holder.
- Stable and light.
- Cushioned seat and flexible backrest.
- Careful finishing: brakes built into the frame, screw cover on the seat...
- Ergonomic handles with adjustable height.
- Parking brakes.
- Curb climber.
- Reflectors on the front wheels and reflective strip on the bag.
- Maximum user weight: 130 kg.
- 2-year guarantee.



### **Technical specifications**



For more technical specifications, please refer to the Rollator New Design manual.



#### **Video**



Availability of the Rollator New Design might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. The medical device mentioned on this document is CE marked according to the European council directive 93/42/EEC and its relatives or the Regulation 2017/745 on medical devices.

The product mentioned in this document is a CE class I device. Please contact Thuasne should you need any additional information on devices classification. Please read carefully the instructions for use, indications and contraindications of the product.

Updated May 2020.

